


Anterior/Posterior Disposable Vitrectomy Cutters

20g - 23g - 25g - 27g

Vitrectomy Cutters



Color codes:


Yellow 20G

Blue 25G

Orange 23G

Violet 27G

Adapters:



AK-0341



AK-0340/7

ALL VIT CUTTERS CAN BE REQUESTED WITH DUAL CUTTING BLADE

INSTRUMENT	Max CPM	20G	23G	25g	27g
GEUDER Megatron S3	2000	AK-0200/BA	AK-0230/BA	AK-0260/BA	AK-0286/BA
ALCON Legacy	800	AK-0201/AA	AK-0231/AA		
AMO Sovereign	800	AK-0202/AA	AK-0232/AA		
AMO Signature	2500	AK-0203/AB	AK-0233/AB		
B&L Stellaris	5000	AK-0205/BF	AK-0235/BF	AK-0265/BF	AK-0287/BF
RUCK Pentasys/2	3000	AK-0206/AA	AK-0236/AA	AK-0266/AA	AK-0288/AA
GEUDER Megatron S4	6000	AK-0207/BA	AK-0237/BA	AK-0267/BA	AK-0289/BA
RUCK Qube	5500	AK-0208/AB	AK-0238/AB	AK-0268/AB	AK-0290/AB
OERTLI	6000	AK-0209/EAA	AK-0239/EAA	AK-0269/EAA	AK-0291/EAA
OPTIKON (Assistant-Pulsar-Antares)	2000/2500	AK-0210/AB	AK-0240/AB	AK-0270/AB	AK-0292/AB
AMO Gemini					
OPTIKON Revolution	6000	AK-0211/AB	AK-0241/AB	AK-0271/AB	AK-0293/AB
ZEISS Visalis	4000/6000				
Millennium	1500	AK-0215/AB	AK-0245/AB	AK-0275/AB	AK-0294/AB
ALCON Accurus	2500	AK-0220/AC	AK-0250/AC	AK-0280/AC	AK-0295/AC
Dorc Associate	2500	AK-0220/AC	AK-0250/AC	AK-0280/AC	AK-0296/AC
Dorc Associate 6000	6000	AK-0219/AC	AK-0249/AC	AK-0279/AC	AK-0297/AC
Dorc Eva	8000	AK-0221/GC	AK-0251/GC	AK-0281/GC	AK-0298/GC
ALCON Constellation*	7500	AK-0225/ABB	AK-0255/ABB	AK-0285/ABB	AK-0299/ABB
ALCON Infiniti	2500	AK-0226/AD	AK-0256/AD		

*needs adapters

Connectors /Aspiration-Pressure line



AK-0398

Pressure line adapter, connector A-D
1 per box (non-sterile)
AVAILABLE DIFFERENT COMBINATIONS

Choose your customized connector for your aspiration line and pressure line

