

Kanülen, etc.

PREMIUM LINIE













MEDIZIN • TECHNIK



Hub Color Code

Hurricane Medical hubs comply with ISO 6009 Color Coding and ISO 594 Conical Fittings with a 6% luer taper to assure use with syringes or devices.

Color		Gauge	ID (in)	ID (mm)	OD (in)	OD (mm)
Natural* (Clear)		30G	.006	.15	.012	.30
Grey		27G	.008	.20	.016	.40
Brown		26G	.010	.25	.018	.45
Orange		25G	.010	.25	.020	.50
Purple		24G	.012	.30	.022	.55
Blue		23G	.013	.33	.025	.60
Black		22G	.016	.41	.028	.70
Green		21G	.020	.51	.032	.80
Yellow		20G	.024	.60	.036	.90
Cream		19G	.028	.70	.043	1.10

*Not an ISO color code

We care!

Kanülen, etc.

PREMIUM LINIE



Anesthesia

For administering local anesthetic prior to ophthalmic surgery
Do not recap needle after use. Discard all needles according to recognized procedures.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Atkinson Tip Anesthesia Needles (single facet with a 22° angle point)				
	1425	Peribulbar Needle - Atkinson 25G x 7/8 in (.50 x 22mm)	10	15,84 €
	1427	Peribulbar Needle - Atkinson 27G x 7/8 in (.40 x 22mm)	10	15,84 €
	1023	Anesthesia Needle - Atkinson 23G x 1 ¼ in (.60 x 32mm) May also be used for peribulbar or retrobulbar administration	10	15,84 €
	1025	Anesthesia Needle - Atkinson 25G x 1 ¼ in (.50 x 32mm) May also be used for peribulbar or retrobulbar administration	10	15,84 €
	1027	Anesthesia Needle - Atkinson 27G x 1 ¼ in (.40 x 32mm) May also be used for peribulbar or retrobulbar administration	10	15,84 €
	1123	Retrobulbar Needle - Atkinson 23G x 1 ½ in (.60 x 38mm)	10	15,84 €
	1125	Retrobulbar Needle - Atkinson 25G x 1 ½ in (.50 x 38mm)	10	15,84 €
	1127	Retrobulbar Needle - Atkinson 27G x 1 ½ in (.40 x 38mm)	10	15,84 €
Sharp Tip Needles (trifacet with a 12° angle point)				
	1225	Retrobulbar Needle - Sharp 25G x 1 ½ in (.50 x 38mm)	10	11,35 €
	1230	Facial Infiltration Needle-Sharp 30G x 5/8 in (.30 x 16mm)	10	11,35 €
Sub-Tenon's Anesthesia (for administering anesthetic agents into the posterior sub-Tenon's space)				
	1320	Sub-Tenon's Cannula 20G x 1 in (.90 x 25mm)	10	50,75 €
	1327	Sub-Tenon's Cannula 27G x 7/8 in (.40 x 22mm)	10	50,75 €
	1330	Sub-Tenon's Cannula 30G x 7/8 in (.30 x 22mm)	10	50,75 €

We care!

Kanülen, etc.

PREMIUM LINIE



Irrigation / Air Injection

For maintaining and forming the anterior chamber by injecting or removing air, fluids, viscoelastics and intraocular medications. Also used to keep the cornea moist. Cannula tip is specially rounded, extra smooth for safe insertion through wound and manipulation in the anterior chamber.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
	2021	Irrigating Cannula - Rycroft 21G x 7/8 in (.80 x 22mm) angled 4mm from end	10	15,84 €
	2023	Irrigating Cannula - Rycroft 23G x 7/8 in (.60 x 22mm) angled 4mm from end	10	15,84 €
	2025	Irrigating Cannula - Rycroft 25G x 7/8 in (.50 x 22mm) angled 4mm from end	10	15,84 €
	2027	Irrigating Cannula - Rycroft 27G x 7/8 in (.40 x 22mm) angled 4mm from end	10	15,84 €
	2030	Irrigating Cannula - Rycroft 30G x 7/8 in (.30 x 22mm) angled 4mm from end	10	15,84 €
	2123	Viscoelastic Cannula 23G x 7/8 in (.60 x 22mm) angled 8mm from end	10	15,84 €
	2125	Viscoelastic Cannula 25G x 7/8 in (.50 x 22mm) angled 8mm from end	10	15,84 €
	2127	Irrigating Cannula 27G x 7/8 in (.40 x 22mm)	10	15,84 €
	2219	Irrigating Cannula - Bishop Harmon 19G x 1 in (1.10 x 25mm) angled from 8mm from end	10	15,84 €
	2220	Irrigating Cannula-Curved 20G x 7/8 in (.90 x 22mm)	10	15,84 €
	2227	Viscoelastic Cannula 27G x 7/8 in (.40 x 22mm) angled 9mm from end	10	15,84 €
	2323	Viscoelastic Cannula 23G x 7/8 in (.60 x 22mm) angled 10mm from end	10	15,84 €
	2327	Viscoelastic Cannula 27G x 7/8 in (.40 x 22mm) angled 10mm from end	10	15,84 €
	2427	Viscoelastic Cannula 27G x 7/8 in (.40 x 22mm) formed 11mm from end	10	19,63 €

We care!

Kanülen, etc.

PREMIUM LINIE



MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
	2520	Anterior Chamber Maintainer 20G x 9/64 in (.90 x 3.5mm) Flexible tubing 10 in (25cm): ID - 1mm OD - 2mm Unique annular rings start .85mm from bevelled tip. Self retaining, does not require suturing	5	134,65 €
	2527	Irrigating Cannula 27G x 7/8 in (.40 x 22mm) angled 7mm from end	10	15,84 €
	2623	Olive Tip Cannula 23G x 7/8 in (.60 x 22mm) curved .33mm end port Olive shape tip is 2mm long x 1mm wide	10	104,98 €
	2727	Irrigating Cannula 27G x 7/8 in (.40 x 22mm) angled 5mm from end	10	15,84 €
	2730	Irrigating Cannula 30G x 7/8 in (.30 x 22mm) angled 5mm from end	10	15,84 €
	2827	Irrigating Cannula 27G x 7/8 in (.40 x 22mm) angled 6mm from end	10	15,84 €



Capsulotomy

Sharp tip cystotome for creating an opening into the anterior capsule, nucleus dislocating or other intraocular applications. Tip angle 90° unless otherwise noted. Each cystotome is packaged to ensure product quality and provide aseptic transfer.

	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
	3023	Irrigating Cystotome 23G x 5/8 in (.60 x 16mm)	10	15,84 €
	3025	Irrigating Cystotome - Pearce 25G x 5/8 in (.50 x 16mm)	10	15,84 €
	3027	Irrigating Cystotome 27G x 5/8 in (.40 x 16mm)	10	15,84 €
	3030	Irrigating Cystotome 30G x 5/8 in (.30 x 16mm)	10	15,84 €

We care!

Kanülen, etc.

PREMIUM LINIE



MEDIZIN • TECHNIK





Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
	3725	Side Cutting Cystotome 25G x 5/8 in (.50 x 16mm)	10	15,84 €
	3925	Irrigating Cystotome 25G x 5/8 in (.50 x 16mm) 45° angled tip	10	19,15 €
Reverse Tip Cystotomes (Tip is bent towards port for downward irrigation)				
	3125	Irrigating Cystotome - Reverse 25G x 5/8 in (.50 x 16mm)	10	17,70 €
	3425	Blunt Tip Cystotome - Reverse 25G x 7/8 in (.50 x 22mm)	10	22,74 €
Short Curve Cystotomes (for capsulotomy on deep-set eyes or small pupils)				
	3225	Irrigating Cystotome 25G x 5/8 in (.50 x 16mm)	10	15,84 €
	3227	Irrigating Cystotome 27G x 5/8 in (.40 x 16mm)	10	15,84 €
Capsulorhexis Cystotomes (Special tip allows control of flap during continuous curvilinear capsulotomies (CCC))				
	3325	Capsulorhexis Cystotome 25G x 5/8 in (.50 x 16mm)	10	26,32 €
	3327	Capsulorhexis Cystotome 27G x 5/8 in (.40 x 16mm)	10	26,32 €
	3627	Capsulorhexis Cystotome - Straight 27G x 5/8 in (.40 x 16mm)	10	26,32 €
Angled and Curved Cystotomes				
	3427	Angled Cystotome 27G x 5/8 in (.40 x 16mm)	10	24,18 €
	3625	Angled Cystotome 25G x 5/8 in (.50 x 16mm)	10	24,18 €
	3825	Endocapsular Cystotome 25G x 5/8 in (.50 x 16mm)	10	16,25 €
	3827	Endocapsular Cystotome 27G x 5/8 in (.40 x 16mm)	10	16,25 €

We care!

Kanülen, etc. PREMIUM LINIE










MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Straight Cystotomes				
	3523	Irrigating Cystotome 23G x 5/8 in (.60 x 16mm)	10	15,84 €
	3525	Irrigating Cystotome 25G x 5/8 in (.50 x 16mm)	10	15,84 €
	3527	Irrigating Cystotome 27G x 5/8 in (.40 x 16mm)	10	15,84 €
	3530	Irrigating Cystotome 30G x 5/8 in (.30 x 16mm)	10	15,84 €



Hydrodissection

For separating the nucleus from the capsule and cortical attachments.
Also for separating and softening layers of nucleus prior to phacoemulsification.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Flat End (flattened for easy placement under anterior capsule and provides a fan-like Distribution of fluid for more effective dissection)				
	4025	Hydrodissector - Pearce 25G x 7/8 in (.50 x 22mm) angled 8mm from end	10	17,77 €
	4027	Hydrodissector 27G x 7/8 in (.40 x 22mm) angled 10mm from end	10	17,77 €
	4030	Hydrodissector 30G x 7/8 in (.30 x 22mm) angled 10mm from end	10	17,77 €
	4125	Hydrodissector - Daya 25G x 7/8 in (.50 x 22mm) formed 11mm and flattened from bend to end	10	17,77 €
	4127	Hydrodissector 27G x 7/8 in (.40 x 22mm) flattened on horizontal plane	10	17,77 €
	4227	Hydrodissector 27G x 7/8 in (.40 x 22mm) flattened on vertical plane	10	17,77 €
	4927	Hydrodissector 27G x 7/8 in (.40 x 22mm) angled 8mm from end	10	17,77 €

We care!

Kanülen, etc.

PREMIUM LINIE



MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
J-shaped (for sub-incisional hydrodissection)				
	4325	J Hydrodissector 25G x 7/8 in (.50 x 22mm) Hook 3mm long x 2.5mm wide	10	40,88 €
	4327	J Hydrodissector 27G x 7/8 in (.40 x 22mm) Hook 3mm long x 2.5mm wide Micro tip for easy insertion to provide sub-incisional hydrodissection	10	40,88 €
	4527	J Hydrodissector - Newman 27G x 7/8 in (.40 x 22mm)	10	40,88 €
	4625	HydroExpression Cannula - Corydon 25G x 7/8 in (.50 x 22mm)	10	62,62 €
	4627	HydroExpression Cannula - Corydon 27G x 7/8 in (.40 x 22mm)	10	62,62 €
Round End (cannula tip is specially rounded, extra smooth for safe insertion through wound and manipulation in the anterior chamber)				
	4427	Hydrodissector 27G x 7/8 in (.40 x 22mm) angled 10mm from end	10	17,56 €
	4430	Hydrodissector 30G x 7/8 in (.30 x 22mm) angled 12mm from end	10	21,70 €
	4725	Hydrodissector 25G x 7/8 in (.50 x 22mm) angled 1mm from end	10	22,74 €
	4727	Hydrodissector 27G x 7/8 in (.40 x 22mm) angled 1mm from end	10	22,74 €
	4827	Hydrodissector - Frankfurt Model 27G x 7/8 in (.40 x 22mm)	10	19,15 €
	5027	Hydrodissector - Chang 27G x 7/8 in (.40 x 22mm) Flat tip angled 90°	10	62,62 €

We care!










Kanülen, etc.

PREMIUM LINIE



Cortex Removal

For removing cortical debris and/or viscoelastic substances.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Bimanual I and A (insert through two side port incisions or side port and phacoemulsification incisions)				
	5121	Bimanual Irrigation Cannula 21G x 7/8 in (.80 x 22mm) .51mm open end irrigation port	10	60,69 €
	5123	Bimanual Aspiration Cannula 23G x 7/8 in (.60 x 22mm) .4mm top port Micro etched 2mm from end for capsul polishing	10	66,21 €
Cortex Extractors (for aspirating when chamber is maintained with a Viscoelastic, anterior chamber maintainer or other independent source)				
	5021	Cortex Extractor 21G x 7/8 in (.80 x 22mm) .4mm top port	10	60,69 €
	5023	Cortex Extractor 23G x 7/8 in (.60 x 22mm) .4mm top port	10	60,69 €
	5323	Cortex Extractor - Charleux 23G x 7/8 in (.60 x 22mm)	10	58,82 €
J Shaped (useful for retracting iris and aspirating residual sub-incisional cortex)				
	5223	J Cannula 23G x 7/8 in (.60 x 22mm) .33mm end port; Hook 2mm long x 2.25mm wide	10	25,36 €
	5225	J Cannula 25G x 7/8 in (.50 x 22mm) .25mm end port; Hook 2mm long x 1.75mm wide; Ideal for small incision surgery	10	25,36 €
	5425	J Cannula - Angled 25G x 7/8 in (.50 x 22mm) Angled 45°, 10mm from end; Hook 2mm long x 1.75mm wide	10	31,43 €
	5425R	J Cortex Extractor - Angled Right 25G x 7/8 in (.50 x 22mm) Angled 45°, 5mm from end J shaped angled right; Hook 1.5mm long x 1.75mm wide	10	35,29 €




We care!

Kanülen, etc.

PREMIUM LINIE




MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Simcoe I and A (Irrigation is provided through cannula with open end Port. Aspiration is provided through 10 inch (25cm) flexible tubing attached to cannula with .4mm top aspiration port)				
	5523	I and A Cannula - Simcoe 23G/23G (.60/.60mm) - Metal Hubs	5	260,92 €
	5623	I and A Cannula - Simcoe 23G/23G (.60/.60mm)	10	201,58 €
	5524	Replacement Tubing 10 in (25cm) ID - 1mm; OD - 2mm; non-sterile	10	17,77 €



Capsule Polishing

For polishing and cleaning residual cortex from the capsular bag in preparation for insertion of the intraocular lens.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Soft Tip (Soft tip for gentle polishing and provides atraumatic insertion and removal)				
	6023	Soft Tip Cannula 23G x 7/8 in (.60 x 22mm) Angled 8mm from end 10mm sleeve	10	39,02 €
	6025	Soft Tip Cannula 25G x 7/8 in (.50 x 22mm) Angled 8mm from end 10mm sleeve	10	39,02 €
	6027	Soft Tip Cannula 27G x 7/8 in (.40 x 22mm) Angled 8mm from end 10mm sleeve	10	39,02 €
	6427	Posterior Capsule Polisher - Engels 27G x 7/8 in (.40 x 22mm); Curved, 23mm sleeve; Beveled down for polishing posterior capsule	10	46,68 €

We care!

Kanülen, etc.

PREMIUM LINIE



MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Micro Etched (Slightly roughened for gentle capsule polishing)				
	6121	Capsule Polisher - Kratz 21G x 7/8 in (.80 x 22mm) .4mm top port; Angled, 9mm from end; micro etched 3mm from end	10	69,45 €
	6123	Capsule Polisher - Kratz 23G x 7/8 in (.60 x 22mm) .4mm top port; Angled, 9mm from end; micro etched 3mm from end	10	69,45 €
	6221	Capsule Polisher 21G x 7/8 in (.80 x 22mm) Angled, 4mm from end; micro etched 5mm from end	10	29,91 €
	6321	Capsule Polisher 21G x 1 ¼ in (.80 x 32mm) Angled, 5mm from end; micro etched 6mm from end	10	29,91 €
	6323	J Capsule Polisher - Thornton 23G x 7/8 in (.60 x 22mm); Hook 2mm long x 2mm wide; Lightly roughened on back of J end for gentle capsule polishing	10	84,15 €
	6423	Olive Tip Polisher 23G x 7/8 in (.60 x 22mm) Angled, 10mm from end; micro etched olive shape tip is 2mm long x 1 mm wide	10	127,27 €

We care!

Kanülen, etc.
PREMIUM LINIE









MEDIZIN • TECHNIK



IOL and Nucleus Manipulation

IOL manipulator for positioning intraocular lenses. Nucleus choppers used during phacoemulsification.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE	
IOL Manipulators (blunt tips avoid damage to intraocular lenses)					
		6801	IOL Manipulator - Lester. 25mm diameter tip Angled, 10mm from tip	5	277,00 €
		6802	IOL Manipulator - Sinsky. 20mm diameter holes and larger Angled, 15mm from tip	5	104,50 €
		6804	IOL Manipulator - Sinsky. 20mm diameter holes and larger	5	104,50 €

We care!


Kanülen, etc.


PREMIUM LINIE



AG

MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Nucleus Choppers (Used to break-up the nucleus within the capsular bag during two-handed phaco techniques)				
	6810	Nucleus Chopper - P. Koch Blunt, round tip; Angled, 6mm from tip; Tip angled 90°, 1.5mm from end	5	125,20 €

	6811	Nucleus Chopper - Beckingsale Blunt, round tip; Angled, 1mm from tip; Tip angled 90°, 1.5mm from end	5	125,20 €
---	------	--	---	----------

We care!

Kanülen, etc.

PREMIUM LINIE








MEDIZIN • TECHNIK



Oculoplastic

For surgical treatment, repair and irrigation of the nasolacrimal system.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Lacrimal Intubation and DCR				
	7001	Lacrimal Intubation Set; Stainless steel probes: 23G x 4 ½ in (.60mm x 11cm); Silicone Tubing Length: 14 in (35cm) ID: .012 in (.30mm) OD: .025 in (.60mm)	5	208,00 €
	7002	Lacrimal Intubation Tubing; Silicone Tubing; Length: 14 in (35cm) ID: .012 in (.30mm) OD: .025 in (.60mm)	10	69,86 €
	7011	DCR Set; Stainless steel probes: 20G x 1 ¾ in (.90 x 45mm) Silicone rod; Length: 15 ¾ in (40cm) Diameter: .032 in (.80mm)	1	37,23 €
	7012	DRC Set - Angled; Angled 45°, 15mm from end; Stainless steel probes: 20G x 1 ¾ in (.90 x 45mm); Silicone rod; Length: 15 ¾ in (40cm); Diameter: .032 in (.80mm)	1	37,23 €
Lacrimal Irrigation				
	7021	Lacrimal Cannula 21G x 1 ¼ in (.80 x 33mm)	10	21,70 €
	7026	Lacrimal Cannula 26G x 1 1/8 in (.45 x 28mm)	10	21,70 €

We care!




Kanülen, etc.

PREMIUM LINIE



Refractive

Unique products designed for LASIK surgery and other types of refractive surgery.













Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Lacrimal Intubation and DCR				
	7121	LASIK Cannula - Daya 21G x 7/8 in (.80 x 22mm) Formed 11mm and flattened from bend to end	10	24,12 €
	7123	LASIK Cannula - Daya 23G x 7/8 in (.60 x 22mm) Formed 11mm and flattened from bend to end	10	24,12 €
	7125	LASIK Irrigator - Rosen 25G x 7/8 in (.50 x 22mm) Formed 9mm and flattened at tip	10	24,12 €
	7127	LASIK Irrigator - 4 Ports 27G x 7/8 in (.40 x 22mm) Formed 9mm from end, tip flattened 1.5mm from end Two .15mm ports 2mm and 4mm from end	10	53,10 €
	7223	LASIK Spatula/Irrigator 23G x 7/8 in (.60 x 22mm) Formed 9mm from end	10	58,82 €
	7225	LASIK Irrigator - 6 Ports 25G x 7/8 in (.50 x 22mm) Formed 7mm from end; Two .15mm ports 2mm, 4mm and 6mm from end	10	53,10 €
	7227	LASIK Irrigator - 4 Ports 27G x 7/8 in (.40 x 22mm) Formed 9mm from end Two .15mm ports 2mm and 4mm from end	10	48,96 €
	7323	LASIK Irrigator - 4 Ports 23G x 7/8 in (.60 x 22mm) Formed 11mm from end Two .2mm ports 2mm and 4mm from end; Flattened tip	10	53,10 €
	7327	LASIK Irrigator - Probst; 27G x 7/8 in (.40 x 22mm) Angled 90°, 7mm from end. Smooth tip, flattened 2mm from end 90° bend allows surgeon to maintain the comfortable superior position while replacing flap	10	24,12 €
	7427	LASIK Irrigator - Formed; 27G x 7/8 in (.40 x 22mm) Formed 11mm from end	10	21,42 €
	7727	LASIK Irrigator - 6 Ports; 27G x 7/8 in (.40 x 22mm) Formed 9mm from end; Two .15mm ports 2mm, 4mm and 6mm from end	10	53,10 €
	7527	LASIK Irrigator - Assis - 3 Ports; 27G x 7/8 in (.40 x 22mm) Formed 9mm from end; Tip flattened 1.5mm from end One .15mm port 2mm, 4mm and 6mm from end	10	53,10 €
	7627	LASIK Spatula / Irrigator; 27G (.40mm); Formed 9mm	10	58,82 €

Kanülen, etc.

PREMIUM LINIE



MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
LASIK Drains (Drains fluid from the orbit and minimizes fluid buildup in corneal bed. Lint-free PVA, soft when moistened)				
	7301	LASIK Drain - Chayet 10mm opening	20	26,12 €
	7311	LASIK Drain - Chayet 12mm opening	20	26,12 €
	7317	LASIK Drain - Gimbel/Chayet 8.5mm opening Unique design to secure flap.	20	32,05 €
Lasik Shields, Spear and Sponges				
	7302	LASIK Shield; Lint free PVA, 8mm diameter with flat edge Protects and moistens corneal flap	20	21,77 €
	7305	LASIK Sponge - Probst; 2 x 4mm, Lint free PVA These convenient small sponges can be used dry or soaked in a vasoconstricting eyedrop to tamponade peripheral neovascular bleeding during LASIK	20 x 5	16,70 €
	7306	LASIK Spear; Lint free PVA, Useful for tissue manipulation and fluid absorption to facilitate flap alignment and adhesion	15 x 5	25,56 €
	7309	LASIK Incise Film; Overall: 4 x 5 in (10 x 12.5cm) Adhesive: 4 x 3 1/2 in (10 x 8.75cm) Mini size perfect for LASIK and other refractive or minor office procedures; Cover lashes, lids and brow	20	37,23 €
Pre-Inked Pads and Pens (Used to coat corneal markers to imprint pattern on the cornea and remains for duration of surgery)				
	7303	Marker Pad - Gentian Violet; Made from lint-free, thick, soft material; Minimizes risk of damaging delicate instruments	10	39,00 €
	7304	Marker Pad - Methylene Blue; Made from lint-free, thick, soft material; Minimizes risk of damaging delicate instruments	10	39,00 €
Marker Pens				
Super moist - will not dry out during procedure Chisel Tip - broad area to mark delicate instruments - fine area for precise lines Fine Tip - for precise lines or to make "pin-point" marks				
	7313	Marker Pen - Gentian Violet	20	39,00 €
	7314	Marker Pen - Methylene Blue	20	39,00 €
	7316	Marker Pen - Gentian Violet, Fine Tip	20	39,00 €

We care!