

Verbrauchsmaterial

PREMIUM LINIE



MEDIZIN • TECHNIK



Vitreoretinal

Products designed specifically for posterior segment surgical procedures.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Infusion (Facilitates posterior segment infusion)				
	7501	Infusion Cannula 2.5mm 20G (.90mm) Flexible Tubing: 10 in (25cm)	10	138,31 €
	7502	Infusion Cannula 4mm 20G (.90mm) Flexible Tubing: 10 in (25cm)	10	138,31 €
	7503	Infusion Cannula 6mm 20G (.90mm) Flexible Tubing: 10 in (25cm)	10	138,31 €
High Viscosity (Designed for injection of silicone oil and other viscous fluids. Special high pressure tubing is attached to thinwall cannula. May also be used as standard infusion cannulas)				
	7504	Viscous Infusion Cannula 4mm 20G (.90mm); Tubing: 4 in (10cm)	10	242,50 €
	7505	Viscous Infusion Cannula 6mm 20G (.90mm); Tubing: 4 in (10cm)	10	242,50 €












Kontakt:

W20 Medizintechnik AG, Bruchstücker 111, D-76661 Philippsburg
Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Infusion (Facilitates posterior segment infusion)				
	7501	Infusion Cannula 2.5mm 20G (.90mm) Flexible Tubing: 10 in (25cm)	10	138,31 €
	7502	Infusion Cannula 4mm 20G (.90mm) Flexible Tubing: 10 in (25cm)	10	138,31 €
	7503	Infusion Cannula 6mm 20G (.90mm) Flexible Tubing: 10 in (25cm)	10	138,31 €
High Viscosity (Designed for injection of silicone oil and other viscous fluids. Special high pressure tubing is attached to thinwall cannula. May also be used as standard infusion cannulas)				
 16 cm				
	7506	Viscous Infusion Cannula 4mm 20G (.90mm); Tubing: 6 in (16cm)	10	242,50 €
	7507	Viscous Infusion Cannula 6mm 20G (.90mm); Tubing: 6 in (16cm)	10	242,50 €
Subretinal Fluid (To aspirate fluid from subretinal space. Soft flexible tips provide atraumatic entry through retinal or macular tears or holes)				
	7720	Subretinal Fluid Cannula - Curved; 20G (.90mm) 6mm flexible tip extension; Overall length with tubing 1 5/16 in (34mm)	10	78,97 €
	7721	Subretinal Fluid Cannula; 20G (.90mm) 6mm flexible tip extension; Overall length with tubing 1 5/16 in (34mm)	10	78,97 €
	7820	Subretinal Fluid Cannula - Curved; 20G (.90mm) 3mm flexible tip extension; Overall length with tubing 1 3/16 in (31mm)	10	78,97 €
	7821	Subretinal Fluid Cannula; 20G (.90mm) 3mm flexible tip extension; Overall length with tubing 1 3/16 in (31mm)	10	78,97 €
	7823	Subretinal Fluid Cannula; 23G (.60mm) 3mm flexible tip extension; Overall length with tubing 1 5/16 in (34mm)	10	115,00 €

Kontakt:

W20 Medizintechnik AG, Bruchstücker 111, D-76661 Philippsburg
Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE



MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
	7725	Subretinal Fluid Cannula; 25G (.50mm) 1.5mm flexible tip extension; Super thinwall inner diameter provides flow; equivalent to 22 gauge (.70mm) cannula; Overall length with tubing 1 1/2 in (31.5mm)	10	83,11
	7825	Subretinal Fluid Cannula; 25G (.50mm) 3mm flexible tip extension; Super thinwall inner diameter provides flow equivalent to 22 gauge (.70mm) cannula Overall length with tubing 1 5/16 in (34mm)	10	83,11 €
Micro Tip Cannulas (designed for injecting or aspirating subretinal fluid. Micro tips for working around and through retinal and macular tears or holes with less chance of engaging the retina)				
	7508	Subretinal Aspirator 27G 20G (.90mm) cannula with 27G (.40mm) tip, angled 2.5mm	10	104,50 €
	7509	Subretinal Injector 30G 20G (.90mm) cannula with 30G (.30mm) tip with gentle curve	10	104,50 €
	7510	Subretinal Injector 30G 20G (.90mm) cannula with 30G (.30mm) tip	10	104,50 €
	7511	Subretinal Injector 30G 20G (.90mm) cannula with 30G (.30mm) tip	10	104,50 €
	7423	Vitreoretinal Cannula 23G x 1 1/8 in (.60 x 28mm) Blunt tip cannula	10	105,00 €
	7512	Subretinal Cannula 32G; 20G (.90mm) cannula with 32G x 3/8 in (.23 x 9mm) tip with gentle curve; Longer tip length is ideal for myopic eyes	10	104,50 €
	7513	Subretinal Cannula 32G 20G (.90mm) cannula with 32G (.23 x 9mm) tip, angled 1mm	10	104,50 €
	7514	Vitreoretinal Cannula 30G 20G (.90mm) cannula with 30G (.30mm) tip	10	104,50 €
	7515	Vitreoretinal Cannula 30G 20G (.90mm) cannula with 30G (.30mm) tip, angled 2.5mm	10	104,50 €
	7516	Vitreoretinal Aspirator 27G 20G (.90mm) cannula with 27G (.40mm) tip	10	104,50 €
	7521	Subretinal Cannula 32G; 20G (.90mm) cannula with 32G x 1/4 in (.23 x 6mm) tip with gentle curve	10	104,50 €

Kontakt:

W20 Medizintechnik AG, Bruchstücker 111, D-76661 Philippsburg

Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE



MEDIZIN • TECHNIK

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
	7522	Vitreoretinal Cannula 32G; 20G (.90mm) cannula with 32G (.23mm) tip	10	104,50 €
	7620	Vitreoretinal Cannula: 20G x 1 1/8 in (.90 x 28mm)	10	49,00 €
	7920	Silicon Oil Cannula; 20G x 5/16 in (.90 x 8mm)	10	48,13 €
	7923	Silicone Oil Cannula; 23G x 5 1/16 in (.60 x 8mm) Blunt tip cannula	10	125,00 €
	7520	Vitreoretinal Brush; 20G (.90mm) 3mm flexible tip extension; for atraumatic brushing of retina	10	201,10 €
	7523	Vitreoretinal Brush; 23G (.60mm) 3mm flexible tip extension	10	221,00 €
	7623	Vitreoretinal Micropick 23G x 1 1/4 in (.90 x 32mm), sharp, angled tip	10	145,00 €
	7625	Vitreoretinal Micropick 25G; 20G (.90mm) cannula with 25G (.50mm) sharp, angled tip; Facilitates removal of premacular, epiretinal and other vitreoretinal membranes	10	78,97 €
	7630	Vitreoretinal Micropick 30G; 20G (.90mm) cannula with 30G (.30mm) sharp, angled tip; Facilitates removal of premacular, epiretinal and other vitreoretinal membranes	10	78,97 €
Viscodissectors (for dissection of vitreoretinal membranes using viscoelastic fluid. A high pressure luer lock tubing set is attached to cannula. Tubing is 8 in (20cm) long: ID - .2 in (.50mm); OD - .06 in (1.5mm) and capable of withstanding injections of viscoelastic fluids)				
	7517	Vitreoretinal viscodissector 90°; 20G (.90mm) cannula with 30G (.30mm) tip, flattened; angled 90°, 1mm from end	10	159,70 €
	7518	Vitreoretinal Viscodissector 70°; 20G (.90mm) cannula with 30G (.30mm) tip, flattened; Angled 70°, 1mm from end	10	159,70 €
	7519	Vitreoretinal Viscodissector 45°; 20G (.90mm) cannula with 30G (.30mm) tip, flattened; Angled 45°, 1mm from end	10	159,70 €

Kontakt:









W20 Medizintechnik AG, Bruchstück 111, D-76661 Philippsburg

Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Fluid Absorption and Administration				
	9006	Surgical Sponge – Gills 2 x 8mm – Lint free PVA, useful for applying or removing fluid in a precise location	20 x 15	52,75 €
	9101	Instrument Wipe 3 x 3 in (7.5 x 7.5cm) Lint free PVA for cleaning delicate instruments. Soft when moistened and will not drag on blades or tips	20	23,43 €
	9301	20 Micron Filter Filters out particles and microbiological contamination .22 micron and larger. Male/female luer lock fittings	50	110,23 €
	9302	Tubing Extension Set Length: 10 in (25cm) ID: .04 in (1mm) OD: .08 in (2mm) Male/female luer lock fittings	20	37,50 €
	9303	Male to Male Connector Overall length 29mm Luer slip fittings connect instrument to tubing; Especially useful for bimanual I and A to connect cannula to I/A tubing	50	27,63 €
	9701	Incise Film Overall size: 5 x 9 in (12.7 x 22.8cm) Adhesive area: 5 x 7 ¼ in (12.7 x 18.4cm) Small incise covers lashes, lids and brow	20	69,45 €
	9802	Eye Fluid Wick Highly absorbent lint free PVA wick; Removes fluid and directs it away from surgical area	20	16,94 €
	9803	Surgical Absorbent Kit: 1 – 9101 Instrument Wipe 1 – 9202 Corneal Light Shield 1 – 9802 Eye Fluid Wick Popular microsurgical accessories combined into one package for convenience and cost savings	20	32,88 €










Kontakt:

W20 Medizintechnik AG, Bruchstücker 111, D-76661 Philippsburg
Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Microsurgical Accesories

A variety of products to complement ophthalmic microsurgery and minor office procedures.

Product Image	Artikel	Bezeichnung	Stück / VPE	Preis pro VPE
Eye Protection (protects the eye during or after surgery or trauma)				
	9401	Clear Eye Shield - non sterile	50	20,53 €
	9402	Clear Eye Shield - sterile Smooth edges need to protective cover; Universal for left or right eye	50	32,40 €
	9202	Corneal Light Shield 8mm diameter - Lint free PVA Keeps cornea moist. Protects retina from operating microscope light rays	20	15,42 €
	9203	Corneal Light Shield 6mm diameter - Lint free PVA	20	15,42 €
Surgical Glides (Useful for placement of intraocular lenses into the posterior chamber. Also protects iris from prolapse or trauma during microsurgical procedures)				
	9201	Surgical Glide; 25mm long x 3mm wide; 1.5mm radius, 6 mil thick	50	38,95 €
	9204	IOL Glide - Sheets; 35mm long x 5mm wide 2.5mm radius, 6 mil thick	50	38,95 €
Fluid Absorption and Administration				
	9001	Eye Spear Fast wicking action for quick absorption	15 x 5	19,63 €
	9003	Absorbent Stick - sterile; 2 5/8 in (6.6cm) overall length 3/16 in (5mm) diameter, 45° angled ends; Soft, non-woven material is rigid enough to withstand gentle pressure. Both ends may be used for absorption	10 x 10	24,25 €
	9004	Absorbent Stick - non sterile; 2 5/8 in (6.6cm) overall length 3/16 in (5mm) diameter, 45° angled ends; Soft, non-woven material is rigid enough to withstand gentle pressure. Both ends may be used for absorption	250	35,91 €

Kontakt:

W20 Medizintechnik AG, Bruchstück 111, D-76661 Philippsburg
Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE



HSI Verbrauchsmaterial Premium Line

Artikelnummer	Artikelbezeichnung	Material	Preis
Universal Vitrektome – cutting rate 600-800 cpm / 20 Gauge			
HS32- 1479B (M-20)	Universal Vitrektom „Megatron“, 20psi	Kunststoff	76,50 €
Universal Vitrectors - cutting rate: 600-800 cpm / 23 Gauge			
HS32- 1517 (M-20)	Universal Vitrektom „Megatron“, 20psi	Kunststoff	119,34 €
Universal Vitrectors cutting rate: 1.200 – 2.500 cpm / 20 Gauge			
HS32- 1479	Universal Vitrektom „Megatron“, 20psi	Kunststoff	119,50 €
Universal Vitrectors cutting rate: 1.200 – 2.500 cpm / 23 Gauge			
HS32- 1517	Universal Vitrektom „Megatron“, 20psi	Kunststoff	129 €

Kontakt:

W20 Medizintechnik AG, Bruchstücker 111, D-76661 Philippsburg
Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE



MEDIZIN • TECHNIK

Alle Lichtfibern erhalten Sie für folgende Maschinen in 20G, 23G und 25G:

Dorc	„Harmony Total TTC“
	„Associate“ (triple-port connector)
Bausch & Lomb	“Storz Premiere”
	“Storz Millenium”
Geuder	“Megatron”
Fritz Ruck	“Pentasys Ocutome”
Alcon	“Grieshaber”

Es kann unter verschiedenen „Ausleuchtungsformen“ gewählt werden.
Unter Angabe der Maschine, Gauge und Ausleuchtungsform besteht die Möglichkeit Ihren Wunsch individuell bei uns zu bestellen.

HSI Verbrauchsmaterial Premium Line

Fiber Optic Endo-Illuminator 20 Gauge

Artikelnummer	Artikelbezeichnung	Ausleuchtung	Preis
HS32-1778	Fiber Optic Endo-Illuminator “Ocutome/Pentasys”	Bullet	58,65 €
HS32-1581	Fiber Optic Endo-Illuminator “Megatron”	Light Pipe	58,65 €
HS32-1792	Fiber Optic Endo-Illuminator “Dorc”	Light Pipe	58,65 €
HS32-1826	Fiber Optic Endo-Illuminator “Megatron”	Total View	58,65 €
HS32-1290	Fiber Optic Endo-Illuminator “Ocutome/Pentasys”	Light Pipe	58,65 €
HS32-1789	Fiber Optic Endo-Illuminator “Synerfit/Synergetics”	Light Pipe	58,65 €

Fiber Optic Endo-Illuminator 23 Gauge

HS32-1788	Fiber Optic Endo-Illuminator “Synerfit/Synergetics”	Light Pipe	64,52 €
HS32-1637	Fiber Optic Endo-Illuminator “Geuder Megatron”	Light Pipe	64,52 €
HS32-1639	Fiber Optic Endo-Illuminator “Ocutome/Pentasys”	Light Pipe	64,52 €
HS32-1675	Fiber Optic Endo-Illuminator “Alcon Accurus”	Light Pipe	64,52 €
HS32-1777	Fiber Optic Endo-Illuminator “Ocutome/Pentasys”	Bullet	64,52 €
HS32-1793	Fiber Optic Endo-Illuminator “Dorc”	Light Pipe	64,52 €
HS32-1801	Fiber Optic Endo-Illuminator “Megatron”	Total View	64,52 €
HS32-1802	Fiber Optic Endo-Illuminator “Ocutome/Pentasys”	Total View	64,52 €
HS32-1849	Fiber Optic Endo-Illuminator “Ocutome/Pentasys”	Cannonball	64,52 €
HS32-1947	Fiber Optic Endo-Illuminator „Dorc“	Shield	64,52 €
HS32-1948	Fiber Optic Endo-Illuminator „Dorc“	Cannonball	64,52 €

Adapter für Fiber Optic

HS32-1623	Adapter “Bausch&Lomb” (adapt to Ocutome)	Alum.	97,75 €
HS32-1680	Adapter “Grieshaber” (adapt to Ocutome)	Alum.	87,98 €

Kontakt:

W20 Medizintechnik AG, Bruchstück 111, D-76661 Philippsburg
Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE



MEDIZIN • TECHNIK

Hochauflösende Vitrektomielinsen aus Silikon zum Einmalgebrauch.
Verfügbar in 5 verschiedenen Ausführungen.

Durch den Stabilisierungsring ist eine Selbstfixierung aller Linsen gewährleistet.
Das Aufnähen der Linse oder eines Fixierendes entfällt.

Linse	Flat (Kilp)	Magnifying	Wide Field	Biconcave	30° Prisma
					
Sichtfeld	36°	30°	48°	25°	33° (offset 30°)
Vergößerung	1,0	1,5	0,5	0,8	1,0
Besonderheit	Standard Linse für Ansicht des HA	Vergrößerte Darstellung der Netzhaut-Pathologie	Darstellung wie Flat (Kilp), nur weiter in die Peripherie	Speziell für Flüssigkeits- Luft-Austausch	Betrachtung der Peripherie in einem bestimmtem Achswinkel
Artikelnummer	VFD10	VMD10	VWFD10	VBCD10	V30PD10
Preis (VPE=10/Box)	279,00 €	279,00 €	279,00 €	279,00 €	279,00 €

Kontakt:

W20 Medizintechnik AG, Bruchstück 111, D-76661 Philippsburg

Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!

Verbrauchsmaterial

PREMIUM LINIE



MEDIZIN • TECHNIK

Staubsaugerkanülen/Infusionsleitungen 23ga Premium Linie



M4469- 23ga Infusionsleitung für Silikonölfüllung, 150mm lang



M4499- 23ga Doppelläufige Kanüle (zum Einbringen von PFCL, ICG, Blue etc)



M4478- 23ga Staubsaugerkanüle „stumpf“



M4459- 23ga Staubsaugerkanüle „Silikonbürste 0,6mm“



M4458- 23ga Staubsaugerkanüle „Soft Tip 2mm“

Besonderheiten:	Alle Kanülen/Infusionsleitungen mit LL- Anschluss
VPE	Jeweils 10 Stk./Box
Artikelnummer	Siehe Produktabbildung
Preis M4469:	249,00 €
Preis M4499:	229,00 €
Preis M4478:	89,00€
Preis M4459:	139,00 €
Preis M4458:	139,00 €

Kontakt:

W20 Medizintechnik AG, Bruchstücker 111, D-76661 Philippsburg

Tel. +49 (0)7256 / 9 25 91 00 | Fax +49 (0)7256 / 9 25 91 10 | info@w2o.de | www.w2o.de

We care!